Website: www.sharpshape.com

Sharp Shape

not just any shape

Emai: sharpshape@comcast.net

Voice Commands and Neutral Casting in AOM TOT jr

Sharp Shape, Copyright ©, All rights reserved, October 2026

1. Words to activate the scanning process

Scanner, **Go**, **Scan**, **Start**, or any combination of these words, like: Scanner Go, Start Scan, or any sentence which contains these words, like: Start the Scanner, Let's Go. False triggering words, like Goal, Scam, etc. and any sentence which contains these words. So try to avoid conversation during a scanning process to avoid false triggering.

2. Words to reset a scan in the middle of a scanning process

Reset, **Retry**, **Redo**, or any combination of these words or any sentence which contains these words, like: Reset Scan, Please Redo, Let's Retry. When these commands are issued, the scanning process will go back to the start of the scanning (frame = 0).

3. Words to stop a scan at the end of a scanning process

Stop, **Enough**, **Finish**, **Exist**, **Quit** or any combination of these words, like: Quit Scanning, Exist Scanning. When these commands are issued, the scanning process will stop and show the 3D contours of the foot on screen.

4. Neutral foot casting methods using the AOMS TOT jr

Some users developed methods to aid neutral foot casting when using the AOMS TOT jr. These methods are used for maintaining subtalar joint neutral position during 3D foot scanning for custom orthotics.

A current active device called 'Foot Positioner' can be found here: https://www.coxorthotics.com/virtualfootimaging.html.

Sulcus Stick was introduced and made by a third party. We are not sure whether it is still active. Info can be found from: https://www.sharpshape.com/aoms_tot/Introduction%20to%20the%20Sulcus%20Stick.pdf.

Without the Positioner or the Stick , some users even used a large rubber band and let the patient to hold the foot in neutral with the rubber band in hand. You can use your imaginations to develop your own method.

If you have any questions, please contact us at sharpshape@comcast.net. - *End of File* -