

12891 Lantana Avenue
Saratoga, CA 95070

 **Sharp Shape**
not just any shape

www.sharpshape.com

Tel & Fax : 408-871-1798

Horizontal Stand Option for Sharp Shape Foot Scanner

Copyright©, Sharp Shape, All rights reserves, December 3, 2012

Some customers are looking for a way to use the Sharp Shape foot scanner horizontally. We have made a design as an option to replace the traditional vertical stand. Please refer to the following photo. All cables should be connected to the back of the scanner. There are four lockable rolling wheels as shown.

Note that our scanner glass cannot bear full weight. The foot can be set flat without applying the body weight on it. You need to make some weight tests before using the scanner horizontally.

If you are only interested in a lower vertical stand or with an angle, we have a PDF document “An Alternative Lower Stand for Sharp Shape Foot Scanner” for you.

If you are interested in this stand as shown, please let us know when you order the scanner. We will replace the traditional vertical stand. The price is the same.



If you have any questions, please contact us at sharpshape@comcast.net. - End of File -